

# TEST KIT (ชุดทดสอบ)

## TEST KIT FOR MILK (ชุดทดสอบผลิตภัณฑ์นม)

### PHOSPHATESMO MI TEST PAPER FOR MILK (กระดาษทดสอบการพาสเจอร์ไรซ์)



#### Phosphatesmo MI

REF 906 12

This test strip allows the specific detection of alkaline phosphatase in milk. In the dairy industry the test is used for the quick and easy quality control of pasteurized milk.

The strip is dipped into the milk and incubated at 36 °C. If the paper remains white, the pasteurization is complete. The test should be stored at 2–8 °C to ensure optimal results.

Type:	test paper
Limit of sensitivity:	0.5 % raw milk in pasteurized milk or 300 U/L alkaline phosphatase in UHT milk
Sufficient for:	50 tests
Color change:	white → yellow

❖ 90612 PHOSPHATESMO MI TEST PAPER FOR MILK (50 TEST) “MN”