


PYKNOMETER (ขนาดความถ่วงจำเพาะ)



Pycnometers. Unadjusted with capacity mark indication at the bottom weighed out up to the 2. decimal place.

Capacity ml	Tolerance ±ml	Thermo- meter °C	Division °C	Unadjusted	Adjusted
Acc. to Gay-Lussac, with ground-in capillary stopper. DIN 12797*					
1	0.006	—	—	3.900.001	3.902.001
2	0.006	—	—	3.900.002	3.902.002
5	0.006	—	—	3.900.005	3.902.005
*10	0.008	—	—	3.900.010	3.902.010
*25	0.010	—	—	3.900.025	3.902.025
*50	0.016	—	—	3.900.050	3.902.050
100	0.024	—	—	3.900.100	3.902.100
With thermometer ST 10/19 and side tube ST 7/16 cap. DIN 12809*					
5	0.006	0-35	0.2	3.920.005	3.922.005
*10	0.008	0-35	0.2	3.920.010	3.922.010
*25	0.010	0-35	0.2	3.920.025	3.922.025
*50	0.016	0-35	0.2	3.920.050	3.922.050
100	0.024	0-35	0.2	3.920.100	3.922.100
Acc. to Reischauer, with cap ST 5/13, DIN 12801*					
5	0.006	—	—	—	3.930.005
10	0.008	—	—	—	3.930.010
*25	0.010	—	—	—	3.930.025
*50	0.016	—	—	—	3.930.050
100	0.024	—	—	—	3.930.100
For dairy products					
20	Cheese, Ice cream	0.2	—	3.935.001	—
50	Cheese	0.5	—	3.935.002	—
20	Ice cream	0.2	—	3.935.003	—
8	Milk	0.1	—	3.935.004	—
10	Milk	0.1	—	3.935.005	—
50	Cream	0.5	—	3.935.006	—